

XEO^{SR}


CUTERA[®]



Tailor your treatment.

XEO^{SR} combines the clinical advantages of laser and light-based technologies with the versatility of specifications to offer the most customizable skin rejuvenation treatments on the market.

XEO^{SR} is designed to give practices the power to adapt to patient demand, and address a wide range of aesthetic concerns.



3D Revitalize Therapy

Revitalize your patients' appearance exclusively with Cutera's xeo SR aesthetic solutions: **limelight, laser genesis and titan**

Every person ages at a different rate and in different ways. Some develop redness, age or sun spots. Others have fine lines or wrinkles. Patients may have one, two or all of these conditions. Cutera's 3D Revitalize Therapy takes a three-tiered approach to revitalizing skin appearance for all of these conditions.



Limelight

is a light-based procedure which addresses problems right at the skin's surface - redness, mottled complexions, freckles and sun damage.



Laser Genesis

is been proven to stimulate new collagen for improvement of fine lines and wrinkles and improve overall skin texture and pore size.



Titan

Infrared heat targeting the deep dermis to tighten the skin.



LimeLight
Benign pigmented lesions (sun/ age spots), vascular lesions and color balance





Photos courtesy of Amy Hammins & Mary Boaglio, Flourish Skin and Laser



Genesis

- Improvement of fine lines and wrinkles
- Improvement of overall skin texture and pore size



Photos courtesy of Richard Fox, MD



Titan
Deep volumetric heating of the dermis, offer revitalize therapy procedure



Photos courtesy of Leonardo Rasi, MD

Epidermal to dermal vascular and pigment imperfections. Treats the broadest range of facial and leg veins including deep reticular vessels

xeo^{SR}								
	Technology Type	Wavelength	Max Energy	Pulse Duration	Repetition Rate	Treatment Spot Sizes	Cooling	
	Limelight®	IPL	1100-520 nm	5-30 J/cm ²	Auto	2 Hz	10 x 30 mm	Contact Sapphire
	Laser Genesis	Laser	1064 nm	300 J/cm ²	0.3-300 ms	10 Hz	3,5,7,10 mm	Gold Plated Copper
	Titan®	IR	1800-1100 nm	5-65 J/cm ²	Auto	---	10 x 30 15 mm	Contact Sapphire

Solely Distributed By:

Dubai
Medica Trading LLC
Business Bay Al Khaleej Al Thejari Street, Vision Tower, 38th Floor
P.O.Box: 31135, Dubai - U.A.E.
Tel: + 971 4 43 70 669
Fax: + 971 4 43 70 668
info@medicagroup.com

Abu Dhabi
Medica Trading LLC
Corniche street, Baynunah Tower 2, 3rd floor, Office #47
P.O.Box: 43385, Abu Dhabi - U.A.E.
Tel: + 971 2 6581 619
Fax: + 971 2 6581 612
info@medicagroup.com

Kingdom of Saudi Arabia
Medica/RIC
Exit 10, King Abdullah Road, Sahhab Tower 1, 5th Floor
P.O.Box: 148, Riyadh 11411 - K.S.A.
Tel: + 966 11 46 54 113
Fax: + 966 11 46 30 135
info@medicagroup.com

Lebanon
Medica Pharmaceuticals s.a.l.
Jal el Dib, Main Road Mallah Center, 11th Floor
P.O.Box: 90601 Jdeideh - Beirut, Lebanon
Tel: + 961 4 710 973
Fax: + 961 4 710 977
info-me@medicagroup.com

Singapore
Medica Singapore Pte Ltd
583 Orchard Road
Suite #14-01, The Forum
P.O.Box: 238884, Singapore
Tel: + 65 67 33 96 89
Fax: + 65 67 33 96 91
info-sg@medicagroup.com