


Total Systemic Health Optimisation and Anti-Aging Technology

Total Xen takes photobiomodulation to the next level, enhancing your treatment offerings by providing a total patient wellness experience encompassing the mind, body and soul. The next level rejuvenation booth utilises evidence-based red and near infrared light to deliver a whole body systemic photobiomodulation experience. With the ability to treat both the face and entire body simultaneously, Total Xen provides an easy service upsell for your clients. Impressive treatment times, ranging from 8-15 minutes allows you to increase and drive additional clients and treatments whilst boosting your business profile.



Red 630nm is primarily absorbed by skin cells. It accelerates cell renewal and repair by increasing collagen production and promoting optimal skin health. It is clinically proven to increase the body's natural form of cellular energy (ATP) which kickstarts cells to regenerate faster.



NIR 850nm penetrates deeper and has even more photobiostimulatory effects. Deep tissue penetration of NIR light stimulates a photochemical reaction within cells that triggers collagen formation and produces systemic healing effects. The process increases blood flow and oxygenation, while improving mitochondrial function, energy production, and cell regeneration.



Clinically Proven Benefits



Skin Health



Muscle Recovery



Body Contouring



Mental Clarity



Reduced Inflammation



Athletic Performance



Sexual Health



Sleep Optimisation



Increased Circulation



Strengthened Immunity

ENGINEERED FOR MAXIMUM EFFICACY AND PERFORMANCE

BENEFITS TO YOUR PRACTICE

- Increase clientele and generate more revenue per patient
- Clinically proven results
- ISO13485:2016 accredited medical device manufacturer
- Australian designed and engineered
- Medically CE marked and TGA-listed
- Advanced training and support
- Results driven = more repeat business
- Easy upsell and promotion of treatment packages
- Mixed modality treatment offerings increase return per client
- Less treatment time required for optimal results
- High patient throughput
- Versatile with customisable protocols
- Proven profitability



Phototherapy has been used in medicine for decades and is backed by thousands of published clinical papers



Supports the body's natural self-healing and regeneration processes to promote a cascade of beneficial effects throughout the body



100% Australian-made photobiomodulation devices with TGA registration and certification of the latest standards of ISO 13485:2016

SPECIFICATIONS

Light Source	Tri-Coloured SMD LED Array
Wavelengths	Red 630 +/- 5nm 660 +/- 5nm NIR 850 +/- 5nm
Class	Medical Class 2A (Risk Group 3)
Power Requirement	240V 50Hz 15A
Power Consumption	3.3KVA
Dimensions	Vertical: 1005 (W) x 1112 (L) x 2020 (H) Horizontal: 1005 (W) x 1937(L)x1112(H)
Weight	280kg



Connect with us

Scan in to learn more or take part in our community



@theaestheticbureau



aestheticbureau.com.au